

UROLOGY

REUSABLE INSTRUMENTS











WORKING ELEMENT MONOPOLAR & BIPOLAR (COMBO)

- Monopolar Resectoscope can be used in complex intrauterine operations without any difficulty
- Biopolar Resectoscope are particularly well suited for infertility patients
- Biopolar resection benefits from generation of plasma to achieve outstanding cutting power
- By using HF current several operations like hysteroscopic endometrium ablation, Myoma resection, adhesiolysis or septum dissection can be performed
- An electrolyte free irrigation fluid is used as distension media

















STONE PUNCH








ELECTRODES

HF- Electrodes, monopolar, single stem (for 24 fr)











CY series Flexible Video Cystoscope

Suction function

One-button control to maintain clear view.

Unique anti-drop design to prevent the suction valve from falling off.



*Patent NO.: ZL202022141662.4

Angle lock design

Lock or release the angle of the distal end for accurate diagnosis and treatment.



*Patent NO.: ZL202022141661.X

Multi-function buttons

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release




Tapered tip

Tapered tip design for easier insertion.



Product P	arameter							
Model	Insertion Tube	Working Channel	Working Length	Bending (Up/Down)	Rotation (Left/Right)	Field of View	Depth of field	Suction function
CY-1356	15.6Fr	6.6Fr	380mm	220°/130°	120°	120°	3 - 100mm	Yes
CY-1355	15.6Fr	6.6Fr	380mm	220°/130°	120°	120°	3-100mm	No





UR series Flexible Video Ureterorenoscope

Secondary bending

Bending U/D 275° Secondary bending to improve operation flexibility



Rotation

Improve operation flexibility and accuracy



*Patent NO.: ZL201621392021

Tapered tip

Tapered tip design for easier insertion



Multi-function buttons

- 1) Zoom in/out
- 2 Capture/record
- ③ Freeze/release



Product Parameter

Model	Insertion	Working	Working	Bending	Rotation	Field	Depth
	Tube	Channel	Length	(Up/Down)	(Left/Right)	of View	of Field
UR-1328	7.5Fr	3.6Fr	670mm	275°/275°	120°	120°	3-50mm

Image Processor



Flexible Video Cystonephroscope



Model	Insertion Tube	Channel	Working Length	Bending (Up/Down)	Field of View	Depth of Field
CY-1356	15.6Fr	6.6Fr	380mm	210°/130°	90°	3-50mm